

W CE MADE IN GERMANY



# World Class Devices by Louise Walsh International

---

With over 30 years of experience, Louise Walsh International delivers the knowledge and expertise to provide exceptional fibroblasting treatments with Plasma Pen and flawless needleless enhancements with HyaPenPro.

Proudly Distributed by  
**Advanced Cosmeceuticals**

 **Plasma Pen™**

 **HyaPenPro™**

**1800 242 011**

# Next-Generation Plasma Devices to Lift, Tighten & Shrink Excess Skin

Plasma Pen™ by Louise Walsh International is one of the most advanced non-invasive, nitrogen plasma devices on the market. German engineered and designed, Plasma Pen delivers next-generation plasma delivery, making it the number one solution for aesthetics and medical professionals worldwide.

**Most common treatments include:**

- Upper & Lower Blepharoplasty
- Lip Lines
- Scar Tissue
- Crow's Feet
- Neck Lift
- Forehead Lines
- Acne Scarring
- Tummy Lift

*All Plasma Pen devices are AC powered with no batteries or power packs to impede the effective and consistent delivery of plasma*

**TGA Registered, CE Certified and FDA Registered**

Ultra-fine, screw-in 0.2mm to 0.4mm nano-probes for absolute precision



Highest plasma delivery intensities and ablative properties in our portfolio



Sleek and ergonomic design

Five power and plasma delivery settings offering maximum versatility and offers precise control

Fingertip controlled, on-demand nitrogen plasma delivery, with an LED display for ease of use

Delivers the fastest spray shading with rapid coverage in our range



Plasma Pen **Platinum**

High plasma delivery intensities and ablative properties



Black brushed aluminium body

Two power and plasma delivery settings

Fingertip controlled, on-demand nitrogen plasma delivery, with an LED display for ease of use



Plasma Pen **Ultra**

## Plasma Pen™ Results

Long lasting and natural results, with short downtimes and fast recovery times.

**Treatment: Upper & Lower Blepharoplasty**



Before Treatment

Immediately After

**Treatment: Crow's Feet**



Before Treatment

Immediately After

**Treatment: Crow's Feet**



Before Treatment

Healed Results

**Treatment: Tummy Lift**



Before Treatment

Healed Results

**Treatment: Neck Lift**



Before Treatment

Healed Results



# Uniform & precise penetration – without a needle

HyaPenPro is a hyaluronic acid filler delivery device designed for uniform, flawless and pain-free volume enhancement. It enables the precise penetration of hyaluronic acid fillers and pharmaceutical grade sculpting and enhancement products into the skin – **all without a needle.**



Before Treatment



After Treatment



Before Treatment



After Treatment

Uses high pressure airflow to create a microscopic hole with a diameter of less than 0.1mm – half the size of the hole left by the smallest injecting needles

Delivers a volume of 0.01ml per actuation

Non-invasive needleless device



Seven pressure settings to target the epidermal and dermal layers of the skin

Predictable airflow delivery mechanism with no social downtime

TGA Registered, CE Certified and FDA Registered

Contact us today to find out more

Proudly Distributed in Australia and New Zealand by **Advanced Cosmeceuticals**

1800 242 011

cs@advancedcosmeceuticals.com.au